

COMMD9 Conjugated Antibody

Catalog No: #C31878

Package Size: #C31878-AF350 100ul #C31878-AF405 100ul #C31878-AF488 100ul #C31878-AF555 100ul #C31878-AF594 100ul #C31878-AF647 100ul #C31878-AF680 100ul #C31878-AF750 100ul #C31878-Biotin 100ul #C31878-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	COMMD9 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human COMMD9
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	COMMD9
Other Names	HSPC166; C11orf55; LINC00610
Accession No.	Swiss-Prot#: Q9Y227NCBI Protein#: BC010892
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966). Modulates Na⁺ transport in epithelial cells by regulation of apical cell surface expression of amiloride-sensitive sodium channel (ENaC) subunits (PubMed:23637203).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.