

HIC1 Conjugated Antibody

Catalog No: #C38632

Package Size: #C38632-AF350 100ul #C38632-AF405 100ul #C38632-AF488 100ul #C38632-AF555 100ul #C38632-AF594 100ul #C38632-AF647 100ul #C38632-AF680 100ul #C38632-AF750 100ul #C38632-Biotin 100ul #C38632-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	HIC1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu, Ms
Immunogen Description	Synthetic peptide of human HIC1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	HIC1
Other Names	hic-1; ZBTB29; ZNF901
Accession No.	Swiss-Prot#: Q6X784NCBI Gene ID: NP_001091672
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

This gene functions as a growth regulatory and tumor repressor gene. Hypermethylation or deletion of the region of this gene have been associated with tumors and the contiguous-gene syndrome, Miller-Dieker syndrome. Alternative splicing of this gene results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.