

ZDHC5 Rabbit pAb

Catalog No: #38025



Package Size: #38025-1 50ul #38025-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZDHC5 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ZDHC5.
Other Names	DHHC5;ZNF375;ZDHC5
Accession No.	Swiss-Prot#:Q9C0B5NCBI Gene ID:25921
Calculated MW	78kDa
Storage	Store at -20oC. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Application Details

WB 1:500 - 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.