

# Arr2 (Phospho-Ser366) Antibody

Catalog No: #SAB345P

Package Size: #SAB345P-1 50ul #SAB345P-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Arr2 (Phospho-Ser366) Antibody  |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates.<br>Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications       | Custom antibody   |
| Species Reactivity | Fruit fly   |
| Immunogen Type     | Peptide-KLH   |
| Target Name        | Arr2  |
| Other Names        | Phosrestin-1; Arrestin-2  |
| Accession No.      | uniprot:P19107  |
| Calculated MW      | 45kDa   |
| Formulation        | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.  |
| Storage            | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.   |

## Application Details

Western blotting: 1:500~1:1000

## Background

Probably plays an important role in the photoreceptor transduction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.