

## Pcli1 Antibody

Catalog No: #SAB324



Package Size: #SAB324-1 50ul #SAB324-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pcli1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Hu Ms
Immunogen Type	Peptide-KLH
Target Name	Pcli1
Other Names	PTB-containing, cubilin and LRP1-interacting protein; PID1
Accession No.	uniprot:Q7Z2X4
Calculated MW	28kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Western blotting: 1:500~1:1000

## Background

Increases proliferation of preadipocytes without affecting adipocytic differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.