

## FLASH Antibody

Catalog No: #SAB265



Package Size: #SAB265-1 50ul #SAB265-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FLASH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Ms
Immunogen Type	Peptide-KLH
Target Name	FLASH
Other Names	FLICE-associated huge protein; CASP8-associated protein 2; CASP8AP2
Accession No.	uniprot:Q9WUF3
Calculated MW	219kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Western blotting: 1:500~1:1000

## Background

Participates in TNF-alpha-induced blockade of glucocorticoid receptor (GR) transactivation at the nuclear receptor coactivator level, upstream and independently of NF-kappa-B. Suppresses both NCOA2- and NCOA3-induced enhancement of GR transactivation (By similarity). Required for histone gene transcription and progression through S phase (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.