

GluRIIB Antibody

Catalog No: #SAB406



Package Size: #SAB406-1 50ul #SAB406-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GluRIIB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Fruit fly
Immunogen Type	Peptide-KLH
Target Name	GluRIIB
Other Names	DGluR-IIBoΩ'½oΩ'½DGluRIIBoΩ'½oΩ'½dGluRIIBoΩ'½oΩ'½dglurIIBoΩ'½oΩ'½GluRIIBoΩ'½oΩ'½Glu-RIIB
Accession No.	uniprot:Q9VMP3
Calculated MW	104kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 70% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.