

## Nlg2 Antibody

Catalog No: #SAB411



Package Size: #SAB411-1 50ul #SAB411-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nlg2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Fruit fly
Immunogen Type	Peptide-KLH
Target Name	Nlg2
Other Names	DNlg2oΩ'½oΩ'½Nlg2-RAoΩ'½oΩ'½Nlg-2oΩ'½oΩ'½nl2oΩ'½oΩ'½Neuroligin
Accession No.	uniprot:Q9NGK5
Calculated MW	137kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 74% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.