

NIBAN2 (Phospho-Tyr593) Antibody

Catalog No: #SAB512P

Package Size: #SAB512P-1 50ul #SAB512P-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NIBAN2 (Phospho-Tyr593) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	NIBAN2
Other Names	Meg-3o Ω ' ω Ω ' ω 2Melanoma invasion by ERKo Ω ' ω Ω ' ω 2 MINERVAo Ω ' ω Ω ' ω 2Niban-like protein 1o Ω ' ω Ω ' ω 2C9orf88o Ω ' ω Ω ' ω 2 FAM129B
Accession No.	uniprot:Q96TA1
Calculated MW	84kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 137% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Background

May play a role in apoptosis suppression. May promote melanoma cell invasion in vitro.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.