

CX3CL1 Rabbit Polyclonal Antibody

Catalog No: #29282

Package Size: #29282-1 50ul #29282-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CX3CL1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB IHC IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	CX3CL1;ABCD-3;C3Xkine;CXC3;CXC3C;NTN;NTT;SCYD1;fractalkine
Accession No.	Swiss-Prot#:P78423NCBI Gene ID:6376
Calculated MW	42kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:200 - 1:2000


IHC 1:50 - 1:200

IF 1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded human brain using CX3CL1 Rabbit pAb (29282) at dilution of 1:150.



Immunofluorescence analysis of NIH-3T3 cells using CX3CL1 Rabbit pAb (29282) at dilution of 1:100 . Blue: DAPI for nuclear staining.



Western blot analysis of various lysates, using CX3CL1 antibody at 1:600 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.