

# ZDHHC18 Rabbit Polyclonal Antibody

Catalog No: #54001

Package Size: #54001-1 50ul #54001-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                            |
|-----------------------|------------------------------------------------------------|
| Product Name          | ZDHHC18 Rabbit Polyclonal Antibody                         |
| Host Species          | Rabbit                                                     |
| Clonality             | Polyclonal                                                 |
| Isotype               | IgG                                                        |
| Purification          | Affinity purification                                      |
| Applications          | WB,IHC                                                     |
| Species Reactivity    | Human,Mouse,Rat                                            |
| Immunogen Description | Recombinant fusion protein of human ZDHHC18 (NP_115659.1). |
| Other Names           | ZDHHC18;DHHC-18;DHHC18                                     |
| Accession No.         | Swiss Prot:Q9NUE0GeneID:84243                              |
| Calculated MW         | 27kDa/42kDa                                                |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.    |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.                |

## Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

## Images



Immunohistochemistry of paraffin-embedded human lung cancer using ZDHHC18 at dilution of 1:100 (40x lens).

## Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.