

NR4A2 Rabbit Polyclonal Antibody

Catalog No: #29244



Package Size: #29244-1 50ul #29244-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	NR4A2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic Peptide of human NR4A2.
Other Names	NR4A2;HZF-3;NOT;NURR1;RNR1;TINUR
Accession No.	Uniprot:P43354GeneID:4929
Calculated MW	59kDa/66kDa
SDS-PAGE MW	67KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using .



Western blot analysis of extracts of Mouse brain, using .

## Background

This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. The encoded protein may act as a transcription factor. Mutations in this gene have been associated with disorders related to dopaminergic dysfunction, including Parkinson disease, schizophrenia, and manic depression. Misregulation of this gene may be associated with rheumatoid arthritis. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.