

ZFPM1 Rabbit Polyclonal Antibody

Catalog No: #54498



Package Size: #54498-1 50ul #54498-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZFPM1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ZFPM1 .
Other Names	FOG;FOG1;ZNF408;ZNF89A;ZC2HC11A;ZFPM1
Accession No.	Swiss Prot:Q8IX07GenelD:161882
Calculated MW	104kDa
SDS-PAGE MW	104kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ZFPM1 at 1:1000 dilution.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.