

GTSF1 Rabbit Polyclonal Antibody

Catalog No: #29241

Package Size: #29241-1 50ul #29241-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GTSF1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human GTSF1.
Other Names	FAM112B;GTSF1
Accession No.	Swiss Prot:Q8WW33GeneID:121355
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

IF 1:50 - 1:200

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.