

## LIPG Rabbit Polyclonal Antibody

Catalog No: #54602

Package Size: #54602-1 50ul #54602-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LIPG Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LIPG (NP_006024.1).
Other Names	LIPG;EDL;EL;PRO719
Accession No.	Uniprot:Q9Y5X9GeneID:9388
Calculated MW	39kDa/56kDa
SDS-PAGE MW	60kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using LIPG antibody.

## Background

The protein encoded by this gene has substantial phospholipase activity and may be involved in lipoprotein metabolism and vascular biology. This protein is designated a member of the TG lipase family by its sequence and characteristic lid region which provides substrate specificity for enzymes of the TG lipase family.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.