

C3AR1 Rabbit Polyclonal Antibody

Catalog No: #29249



Package Size: #29249-1 50ul #29249-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	C3AR1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human C3AR1 (NP_004045.1).
Other Names	C3AR1;AZ3B;C3AR;HNFAG09
Accession No.	Swiss Prot:Q16581GeneID:719
Calculated MW	53kDa
SDS-PAGE MW	55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using C3AR1 at 1:3000 dilution.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.