

Recombinant Human Transcription factor GATA-6(GATA6)

Catalog No: #AP0106

Package Size: #AP0106-1 20ug #AP0106-2 100ug #AP0106-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Transcription factor GATA-6(GATA6)
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	Gata binding factor 6; Gata binding protein 6; GATA-binding factor 6; Gata6; GATA6_HUMAN; Transcription factor Gata 6; Transcription factor GATA-6
Tag Info	N-terminal His-tagged
Target Sequence	1-595aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20oC/-80oC. The shelf life of lyophilized form is 12 months at -20oC/-80oC.</p>

Application Details

Greater than 85% as determined by SDS-PAGE.

Background

Transcriptional activator

Note: This product is for in vitro research use only and is not intended for use in humans or animals.