

# Recombinant Human Cation channel sperm-associated protein 1(CATSPER1)

Catalog No: #AP0111

Package Size: #AP0111-1 20ug #AP0111-2 100ug #AP0111-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cation channel sperm-associated protein 1(CATSPER1)
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Sequence	1-445aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

## Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.