

Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit

Catalog No: #EK21006



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	VASHL, Vasohibin-2, Vasohibin-like protein
Accession No.	Uniprot:Q86V25
Storage	6 months

Application Details

Detection Range: 0.78-50ng/mL Assay Type: Sandwich Sensitivity: 0.24ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.