

Human RGS2 / Regulator of G-protein signaling 2 ELISA Kit

Catalog No: #EK21007



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human RGS2 / Regulator of G-protein signaling 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	RGS2, Cell growth-inhibiting gene 31 protein, G0/G1 switch regulatory protein 8, G0S8, GIG31
Accession No.	Uniprot:P41220
Storage	6 months

Application Details

Detection Range: 0.625-40 ng/mL Assay Type: Sandwich Sensitivity: 0.317 ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.