

Human CCDC25 / Coiled-coil domain-containing protein 25 ELISA Kit

Catalog No: #EK21008



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human CCDC25 / Coiled-coil domain-containing protein 25 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:Q86WR0
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.086ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.