

# Human GGCX / Vitamin K-dependent gamma-carboxylase ELISA Kit

Catalog No: #EK21009



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human GGCX / Vitamin K-dependent gamma-carboxylase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Gamma-glutamyl carboxylase,Peptidyl-glutamate 4-carboxylase,Vitamin K gamma glutamyl carboxylase,GC
Accession No.	Uniprot:P38435
Storage	6 months

## Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39.9pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.