

Human MT-CYB / Cytochrome b ELISA Kit

Catalog No: #EK21010

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human MT-CYB / Cytochrome b ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Complex III subunit 3,Complex III subunit III,Cytochrome b-c1 complex subunit 3,Ubiquinol-cytochrome-c reductase complex cytochrome b subunit,COB,CYTB,MTCYB
Accession No.	Uniprot:P00156
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:38pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.