

Mouse Chrac1 / Chromatin accessibility complex protein 1 ELISA Kit

Catalog No: #EK21047



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Chrac1 / Chromatin accessibility complex protein 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	CHRAC-1,DNA polymerase epsilon subunit p15,NF-YC-like protein,YC-like protein 1,YCL1
Accession No.	Uniprot:Q9JKP8
Storage	6 months

Application Details

Detection Range:7.8-500pg/mL Assay Type:Sandwich Sensitivity:3.82pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.