

# Mouse Apela / Apelin receptor early endogenous ligand ELISA Kit

Catalog No: #EK21054

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Apela / Apelin receptor early endogenous ligand ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Protein Elabela, Protein Toddler, Ela, Ende, Gm10664, Tdl, ELA
Accession No.	Uniprot:P0DMC4
Storage	6 months

## Application Details

Detection Range: 0.78-50ng/mL Assay Type: Competitive Sensitivity: 0.397ng/mL

## Product Description

Detection Method: Competitive

Note: This product is for in vitro research use only and is not intended for use in humans or animals.