

PAQR3 (Phospho-Ser8) Antibody

Catalog No: #SAB524P

Package Size: #SAB524P-1 50ul #SAB524P-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAQR3 (Phospho-Ser8) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	PAQR3
Other Names	Progestin and adipoQ receptor family member III,Raf kinase trapping to Golgi,RKTG
Accession No.	uniprot:Q6TCH7
Calculated MW	36kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 142% glycerol.
Storage	Store at -20oΩ½oΩ½C for long term preservation (recommended). Store at 4oΩ½oΩ½C for short term use.

Application Details

Western blotting: 1:500~1:1000

Background

Functions as a spatial regulator of RAF1 kinase by sequestering it to the Golgi.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.