

Recombinant 2019-nCoV S Protein RBD (Mammalian, C-His)

Catalog No: #AG0068

Package Size: #AG0068-1 100ug #AG0068-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant 2019-nCoV S Protein RBD (Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	2019-nCov RBD Protein, 2019-nCoV Spike RBD Protein
Calculated MW	28.03kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.