

# Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)

Catalog No: #AG0071

Package Size: #AG0071-1 100ug #AG0071-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, South African (ZA) variant (B.1.351), Beta
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

## Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.