

# Mouse Sclerostin (SOST) ELISA Kit

Catalog No: #EK20161



Package Size: #EK20161-1 48T #EK20161-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Sclerostin (SOST) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:6.5pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.