

Mouse Osteocalcin (OC) ELISA Kit

Catalog No: #EK20179



Package Size: #EK20179-1 48T #EK20179-2 96t

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Osteocalcin (OC) ELISA Kit
Applications	ELISA
Species Reactivity	Ms

Application Details

Detect Range:62.5-4000pg/ml

Sensitivity:23.58pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.