

Recombinant Severe acute respiratory syndrome coronavirus 2 RNA-dependent RNA polymerase(NSP12)

Catalog No: #AP89627

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89627-1 20ug #AP89627-2 100ug #AP89627-3 1mg

Description

Product Name	Recombinant Severe acute respiratory syndrome coronavirus 2 RNA-dependent RNA polymerase(NSP12)
Brief Description	Recombinant Protein
Host Species	E.coli
Immunogen Description	Expression Region:1-932aaSequence Info:Full Length of Mature Protein
Calculated MW	108.2 kDa
Tag Info	N-terminal His-tagged
Target Sequence	SADAQSFLNRVCGVSAARLTPCGTGTSTDVVYRAFDIYNDKVAGFAKFLKTNCCRFQEKDEDDNLIDSYFVVK RHTFSNYQHEETIYNLLKDCPAVAKHDFKFRIDGDMVPHISRQRLTKYTMADLVYALRHFDEGNCDTLKEILV TYNCCDDDYFNKKDWYDFVENPDILRVYANLGERVRQALLKTVQFCDAMRNAGIVGVLTDNQLNGNWDYD FGDFIQTPGSGVPVVDSSYLLMPILTLTRALTAESHVDTDLTKPYIKWDLKDYFTEERLKLFDYFYKYWDQ TYHPNCVNCLDDRCILHCANFNVLSTVFLPTSFGPLVRKIFVDGVFPVSTGYHFRELGVVHNQDVLNLSRL SFKELLVYAADPAMHAASGNLLDKRTTCFSVAALTNNVAFQTVKPGNFNKDFYDFAVSKGFFKEGSSVELKH FFFAQDGNAAISDYDYRYNLPMTCDIRQLLFVVEVDKYFDCYDGGCINANQVIVNNLDKSAGFPFNKWGKA RLYYDSMSYEDQDALFAYTKRNVIPITITQMNLYAISAKNRARTVAGVSICSTMTNRQFHQKLLKSIAATRGAT VVGTSKFYGGWHNMLKTVYSDVENPHLMGWDPKCDRAMPNMLRIMASLVLRKHHTTCCSLSHRFYRLAN ECAQVLSEMVCMCGSLYVKPGGTSSGDATTAYANSVFNICQAVTANVNALLSTDGNIADKYVRNLQHRLEYE CLYRNRDVTDFVNEFYAYLRKHFSMMILSDDAVVCFNSTYASQGLVASIKNFKSVLYYQNNVFMSEAKCWT ETDLTKGPHEFCSQHTMLVKQGDDYVYLPYPDPSPRILGAGCFVDDIVKTDGTLMIERFVSLAIDAYPLTKHPNQ EYADVFLYLYQYIRKLHDELTHMLDMYSVMLTNDNTSRYWEPEFYEAMYTPHTVLQ
Formulation	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.