

# MX1 / MX2 Polyclonal Antibody Gold Conjugated

Catalog No: #C90405Gold

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MX1 / MX2 Polyclonal Antibody Gold Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC
Species Reactivity	Hu Ms Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human IFI78/MX1
Conjugates	Gold
Target Name	MX1 / MX2
Other Names	MX1 / MX2
Accession No.	Uniprot:P20591Gene ID:4599
Formulation	20 mM TBS (pH=8.0) with 1% BSA and 0.03% Proclin300.
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.