

# PRRSV-GP5 Polyclonal Antibody Bio Conjugated

Catalog No: #C92170Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PRRSV-GP5 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	PRRSV
Immunogen Description	KLH conjugated synthetic peptide derived from PRRSV GP5: 101-200/200
Conjugates	Bio
Target Name	PRRSV-GP5
Other Names	PRRSV-GP5
Accession No.	Uniprot:O55492
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.