

Lenalidomide (CC-5013)

Catalog No: #S1642

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lenalidomide (CC-5013)
Brief Description	Inhibitors
Purification	98.00%
Target Name	TNF inhibitor
Calculated MW	259.26
Formulation	C13H13N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Cancer

SMILES: NC1=CC=CC2=C1CN(C1CCC(=O)NC1=O)C2=O

Pathways: Apoptosis

Receptor: TNF- α

Boiling pt:

Melting pt: 269-271°C

Solubility: DMSO: 35 mg/mL

Appearance: White Powder

Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Muller GW, et al. Bioorg Med Chem Lett, 1999, 9(11), 1625-1630.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.