

## 2,4'-DIHYDROXYBENZOPHENONE

Catalog No: #S7036

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	2,4'-DIHYDROXYBENZOPHENONE
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	214.22
Formulation	C13H10O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:OthersSMILES:C1=CC=C(C(=C1)C(=O)C2=CC=C(C=C2)O)OPathways:OthersReceptor:OthersBoiling pt:Melting  
pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a  
while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.