

WDR43 Rabbit Polyclonal Antibody

Catalog No: #53194

Package Size: #53194-1 50ul #53194-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	WDR43 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of mouse WDR43 (NP_783570.1).
Other Names	WDR43;NET12;UTP5
Accession No.	Uniprot:/GeneID:72515
Calculated MW	75kDa
SDS-PAGE MW	75KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using WDR43 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.