

Recombinant Human Transient receptor potential cation channel subfamily V member 6(TRPV6),partial

Catalog No: #AP89635



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP89635-1 20ug #AP89635-2 100ug #AP89635-3 1mg

Description

Product Name	Recombinant Human Transient receptor potential cation channel subfamily V member 6(TRPV6),partial
Brief Description	Recombinant Protein
Host Species	E.coli
Target Name	His-tagged
Other Names	TRPV6; ECAC2; Transient receptor potential cation channel subfamily V member 6; TrpV6; CaT-like; CaT-L; Calcium transport protein 1; CaT1; Epithelial calcium channel 2
Accession No.	Uniprot No.øO Q9H1D0
Tag Info	N-terminal His-tagged
Target Sequence	Partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.