

Recombinant Human B-cell CLL/lymphoma 7 protein family member A(BCL7A)

Catalog No: #AP89658

Package Size: #AP89658-1 20ug #AP89658-2 100ug #AP89658-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human B-cell CLL/lymphoma 7 protein family member A(BCL7A)
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	B cell CLL/lymphoma 7; B cell CLL/lymphoma 7 protein family member A; B cell CLL/lymphoma 7A; B-cell CLL/lymphoma 7 protein family member A; BCL 7; BCL7; bcl7a; BCL7A_HUMAN
Accession No.	Q4VC05
Target Species	Hu
Tag Info	N-terminal His-tagged
Target Sequence	1-210 full length protein MSGRSVRAET RSRAKDDIKR VMAAIEKVRK WEKKWVTVDG TSLRIYKWVP VTEPKVDDKN KNKKKGKDEK CGSEVTPEN SSSPGMMMDMH DDNSNQSSIA DASPIKQENS SNSSPAPEPN SAVPSDGTEA KVDEAQADGK EHPGAEDASD EQNSQSSMEH SMNSSEKVDR QPSGDSGLAA ETSAISQDLE GVPPSKMKML EASQQNSEEM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.