

Recombinant Human NACHT, LRR and PYD domains-containing protein 3(NLRP3),partial

Catalog No: #AP89639

Package Size: #AP89639-1 20ug #AP89639-2 100ug #AP89639-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NACHT, LRR and PYD domains-containing protein 3(NLRP3),partial
Host Species	E.coli
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Expression Region: 733-1036aa Sequence Info:partial
Target Name	NLRP3
Other Names	(Angiotensin/vasopressin receptor AII/AVP-like)(Caterpillar protein 1.1)(CLR1.1)(Cold-induced autoinflammatory syndrome 1 protein)(Cryopyrin)(PYRIN-containing APAF1-like protein 1)
Accession No.	Q96P20
Target Species	Homo sapiens (Human)
Calculated MW	40.7 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	LFSVLSTSQSLTELDLSDNSLGDPMRVLCETLQHHPGCNIRRLWLGRGCLSHECCFDISLVSSNQKLVELDLSDNALGDFGIRLLCVGLKHLLCNLKKLWLVSCCLTSACQDLASVLSTSHSLTRLVYGENALGDSGVAILCEKAKNPQC�NLQKLGLVNSGLTSVCCSALSSVLSTNQNLTHLYLRGNTLGDKGIKLLCEGLLHPDCKLQVLELDNCNLTSHCCWDLSTLLTSSQSLRKLSLGNNDLGDLDGVMMFCEVLKQQSCLLQNLGLSEMYFNYETKSALETLQEEKPELTVVFEPWS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.