

GLP1R Rabbit Polyclonal Conjugated Antibody

Catalog No: #C29668



Package Size: #C29668-AF350 100ul #C29668-AF405 100ul #C29668-AF488 100ul #C29668-AF555 100ul #C29668-AF594 100ul
#C29668-AF647 100ul #C29668-AF680 100ul #C29668-AF750 100ul #C29668-Biotin 100ul #C29668-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GLP1R Rabbit Polyclonal Conjugated Antibody
Clone No.	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GLP1R (NP_002053.3).
Other Names	GLP1R;GLP-1;GLP-1-R;GLP-1R
Accession No.	Swiss Prot:P43220GeneID:2740
Calculated MW	53kDa
SDS-PAGE MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.