

Monkey Angiopoietin like protein 3 (ANGPTL3) ELISA Kit

Catalog No: #EK20373

Package Size: #EK20373-1 48T #EK20373-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monkey Angiopoietin like protein 3 (ANGPTL3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:1.56-100 ng/mL

Sensitivity:0.53ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.