

## Donkey Anti-Goat IgG (H&amp;L) DyLight 488

Catalog No: #L35657

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Donkey Anti-Goat IgG (H&L) DyLight 488
Brief Description	Donkey
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	Goat
Immunogen Description	Goat IgG
Conjugates	DyLight 488
Formulation	Liquid in 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for six months. Avoid freeze/thaw cycles.

## Application Details

IF (1/100 - 1/800), FC (1/100 - 1/800)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.