

Recombinant Human Tau protein, Tau(2N4R)

Catalog No: #AP90002

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Tau protein, Tau(2N4R)
Brief Description	Recombinant Proteins
Host Species	Mammalian Cell
Purification	Greater than 95% as determined by SDS-PAGE.
Immunogen Description	1-441aa
Target Name	His-tagged
Other Names	Tau-4, 2N4R
Accession No.	P10636-8
Target Species	Human
SDS-PAGE MW	46 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	MAEPRQEFEVMDHAGTYGLGDRKDQGGYTMHQDQEGD TDAGLKESPLQTP TEDGSEEPGSETSDAKSTP TAEDVTAPLVDEGAPGKQAAAQPHT EIP EGT TAAEEAGIGDTPSLEDEAAGHV TQARMVSKSKDGTGSDDKKA KGADGKTKIATPRGAAPP GQKGQANATRI PAKTP PAPKTPPSSGEPPKSGDRSGYSSPGSPGTPGSRSRTPS LTPPTREP KKVAVVRTPPKSPSSAKSRLQTAPVPM PDLKNV KSKIGSTENLKHQPGGGKVQIINKKLDLSNVQ SKCGSKDNIKHVPGGGSVQIVYK PVDLSKVTSKCGSLGNIHHKPGGGQVEVKSEKLD FDKDRVQSKIGSLDNIT HVPGGGNKKIETHKLT FREN AKAKTDHGAEIVYKSPVVS GDTSPRHLSNVSS TGSIDMV DSPQLATL ADEVSA SLAKQGL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.