

β-Actin Rabbit mAb

Catalog No: #52901



Package Size: #52901-1 50ul #52901-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	β-Actin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	beta Actin Antibody detects endogenous levels of total beta-Actin protein.
Immunogen Description	A synthesized peptide derived from human beta Actin
Other Names	ACTIN; ACTB; BETA-ACTIN; beta actin; actin beta;
Accession No.	Swiss Prot:P60709
Calculated MW	42kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:5000~1:20000  
IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of extracts from 3T3,293, HeLa ,Mouse Brain and Rat Brain using beta Actin Rabbit mab Antibody #52901.

Published Papers

el at., TSC22D3 as an immune-related prognostic biomarker for acute myeloid leukemiaInScienceOn2023 Jul 22byYang Li 1 2, Hanying Huang et al..PMID: 37575189, , (2023)  
[PMID:37575189](#)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.