

GAPDH Rabbit mAb

Catalog No: #52902

Package Size: #52902-1 50ul #52902-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GAPDH Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	H, M, R, Mk, Ch, Ze, Fish
Specificity	GAPDH Antibody detects endogenous levels of total GAPDH
Immunogen Description	A synthesized peptide derived from human GAPDH
Other Names	aging-associated gene 9 protein; G3P; G3PD; GAPD; GAPDH; glyceraldehyde 3-phosphate dehydrogenase; Glyceraldehyde-3-phosphate dehydrogenase; MGC88685;
Accession No.	Swiss Prot:P04406
Calculated MW	36kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:5000~1:30000

IHC 1:100~1:500 ICC/IF 1:100~1:250

IP 1:50 FC 1:50

Images



Western blot analysis of extracts from 3T3,293, HeLa ,Mouse Brain and Rat Brain using GAPDH Rabbit mab #52902.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.