

## HARS Rabbit mAb

Catalog No: #57028

Package Size: #57028-1 50ul #57028-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HARS Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HARS Antibody detects endogenous levels of total HARS
Immunogen Description	A synthesized peptide derived from human HARS
Other Names	HARS; HisRS; histidine transylase; HRS; USH3B;
Accession No.	Uniprot:P12081
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of HARS expression in 293T cell lysate.

## Product Description

Plays a role in axon guidance.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.