

# Rabbit Anti-Human IgM

Catalog No: #L35658

Package Size: #L35658-1 100ul #L35658-2 500ul #L35658-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Rabbit Anti-Human IgM   |
| Brief Description     | Rabbit Anti-Human IgM   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications          | E, WB, IH   |
| Species Reactivity    | Human   |
| Immunogen Description | Human IgM   |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide  |
| Storage               | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                  |

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.