

Calreticulin Antibody PE-Cy7 Conjugated

Catalog No: #C06095PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calreticulin Antibody PE-Cy7 Conjugated
Brief Description	Calreticulin Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	FC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Calreticulin
Conjugates	PE-Cy7
Target Name	Calreticulin
Other Names	RO; CRT; SSA; cC1qR; HEL-S-99n; Calreticulin; CRP55; Calregulin; Endoplasmic reticulum resident protein 6; ERp6; HACBP; grp6; CALR; CRTC
Accession No.	#P27797 NCBI Gene ID811
Calculated MW	44kD
Concentration	1mg /ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Calcium-binding chaperone that promotes folding, oligomeric assembly and quality control in the endoplasmic reticulum (ER) via the calreticulin calnexin cycle. This lectin interacts transiently with almost all of the monoglucosylated glycoproteins that are synthesized in the ER. Interacts with the DNA-binding domain of NR3C1 and mediates its nuclear export. Involved in maternal gene expression regulation. May participate in oocyte maturation via the regulation of calcium homeostasis (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.