

Thrombospondin 4 Polyclonal Antibody

Catalog No: #44224



Package Size: #44224-1 50ul #44224-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Thrombospondin 4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC-p;IF/ICC;ELISA
Species Reactivity	Human; Mouse;Rat
Specificity	Thrombospondin 4 Polyclonal Antibody detects endogenous levels of Thrombospondin 4 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human THBS4.
Accession No.	Swiss Prot:P35443GenelD:7060
Calculated MW	106kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000;
IHC-p 1:100-500;IF/ICC 1:100-500;
ELISA 1:5000-20000

Images



Western blot analysis of 293T lysis using Thrombospondin 4 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using THBS4 Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.