

## HRPT2 / CDC73 Rabbit mAb

Catalog No: #57154



Package Size: #57154-1 50ul #57154-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HRPT2 / CDC73 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	HRPT2 / CDC73 Antibody detects endogenous levels of total HRPT2 / CDC73
Immunogen Description	A synthesized peptide derived from human HRPT2 / CDC73
Other Names	Cdc73; FIHP; HPTJT; HRPT1; HRPT2; Hyperparathyroidism 2; HYX;
Accession No.	Q6P1J9
Calculated MW	61kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of HRPT2/CDC73 expression in 293T cell lysate.

## Background

Tumor suppressor probably involved in transcriptional and post-transcriptional control pathways. May be involved in cell cycle progression through the regulation of cyclin D1/PRAD1 expression.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.